



UNDERIVITIZED AMINO ACIDS

Application note # QAN-0036

PARAMETERS


NQAD
GAS: Air
FILTER: 5
GAIN: 10x
EVAPORATION TEMP: 35°C

HPLC
MOBILE PHASE:

Time	5 mM HFBA, pH 1.0 adjusted with 0.1% TFA in Water	ACN
0	100	0
6	100	0
8	85	15
25	35	65

FLOW RATE: 1.2 mL/min
COLUMN: Alltima C18, 5u, 250 x 4.6
mm
COLUMN TEMP: Ambient
INJECTION VOLUME :10 µL

CHROMATOGRAM



1 = Gly, 2 = Ala, 3 = Cys, 4 = Pro, 5 = Val, 6 = Tyr, 7 = Leu, 8 = Trp



**QUANT TECHNOLOGIES,
LLC**
1463 94TH Lane NE
Blaine, MN, USA 55449
763-398-0508
www.Quant-NQAD.com