



VITAMIN D2

Application note # QAN-0042

PARAMETERS

NQAD

GAS: Air
FILTER: 5 sec
GAIN: 2x-1kx
EVAPORATION TEMP: off

HPLC

MOBILE PHASE: 25:75 MeOH:ACN
FLOW RATE: 1.2 mL/min
COLUMN: Alltech Alltima C18, 3u,
4.6 x 33 mm
COLUMN TEMP: Ambient
INJECTION VOLUME :10 µL
CONCENTRATION: 20-20,000 ng on
column
ANALYTE: Vitamin D2



QUANT TECHNOLOGIES, LLC

1463 94TH Lane NE
Blaine, MN, USA 55449
763-398-0508
www.Quant-NQAD.com

CHROMATOGRAMS/DATA

