



# CARBARYL REPRODUCIBILITY

Application note # QAN-0053

## PARAMETERS

### NQAD

GAS: Air  
FILTER: 0.6 sec  
GAIN: 10x  
EVAPORATION TEMP: 35C

### HPLC

MOBILE PHASE: 50:50  
MeOH:Water  
FLOW RATE: 0.5 mL/min  
COLUMN: Vision HT C18 B, 1.5  $\mu$ ,  
50 x 2 mm  
INJECTION VOLUME : 1  $\mu$ L  
ANALYTE: 100 ng oxamyl, 200 ng  
carbaryl



**QUANT TECHNOLOGIES**  
1463 94<sup>TH</sup> Lane NE  
Blaine, MN, USA 55449  
763-398-0508  
[www.Quant-NQAD.com](http://www.Quant-NQAD.com)

